


Selecting Technology Levels Using Inspire Software

Counseling Tools in Inspire software has a Lifestyle Solutions component that can be used to enhance the process of selecting the right technology for patients. It includes a patient questionnaire that helps patients understand the features and benefits of each technology level, especially as it relates to patient-specific needs. Lifestyle Solutions can be accessed through Counseling Tools, from the icon on the Get Started screen or the left navigation bar.



Lifestyle Solutions Icon




Left Navigation Bar

Step-by-Step Process: Lifestyle Solutions


- 1) Instruct the patient to select a rating
- 2) (Rarely/Never, Occasionally, Often, Very Often) for each of the ten different listening environments on the Lifestyle Solutions Questionnaire, based upon their desire to participate in the environment activity with improved hearing.
- 3) Click either the Show Recommendations (show most appropriate products) or Show All (show all products) button to navigate to the Product Family screen.
- 4) Click either Mobility or Focus to view the product tier details on the Product Selection screen. Click on any or all desired products to include in recommendations (if not already selected).



Lifestyle Solutions Questionnaire




Product Family Screen



Product Selection Screen

NOTE: Selected products will be included on the **Summary** screen. Inspire allows the professional to toggle between the wired and wireless product families through a hyperlink on the top of the page.

- 5) Click the **Show Styles** button to display the different hearing instrument styles. To select a particular style to include on the **Summary** screen, click the style picture. A checkmark will appear in the space next to the picture.
- 6) Click the **Show Summary** button to display all of the compiled information. A description of the recommended technology levels, styles and features will be displayed.
- 7) Type customized patient information in the **Notes** box. This will print on the **Patient Report**.
- 8) Click on the **Print Report** button to print the report for the patient.



Styles Screen



Summary Screen