

PERFORMANCE

**MAXIMUM (12dB)**

Channels Strategic Feedback Control SpeechTrax<sup>2</sup> Dynamic Directional Detection

FrequencyTrax Noise Speech Quiet Adaptive Environmental Sequencing

Speech In Noise Machine Noise Wind

**PREMIUM**

**MODERATE (8dB)**

Channels Strategic Feedback Control SpeechTrax<sup>2</sup> Dynamic Directional Detection

FrequencyTrax Noise Speech Quiet Adaptive Environmental Sequencing

Speech In Noise Machine Noise Wind

**ADVANCED**

**LIGHT (6dB)**

Channels Strategic Feedback Control SpeechTrax<sup>2</sup> Dynamic Directional Detection

FrequencyTrax Noise Speech Quiet Adaptive Environmental Sequencing

Machine Noise Wind

**SELECT**

WIRELESS

SurfLink™ Media SurfLink Remote **ST Boost** SurfLink Programmer Synchronized User Controls

Machine Noise Wind **PREMIUM** Binaural Spatial Mapping Dynamic Directional Detection Binaural Telephone Mode

SurfLink™ Media SurfLink Remote **ST Boost** SurfLink Programmer Synchronized User Controls

Dynamic Directional Detection **ADVANCED** Binaural Spatial Mapping Binaural Telephone Mode

SurfLink™ Media SurfLink Remote **ST Boost** SurfLink Programmer Synchronized User Controls

Binaural Telephone Mode

FEATURES

Self Learning Swap Fit Auto Path Voice Indicators

Music Memory Memories Television Memory MT Connect™

Live Real Ear Measurement Self Check T<sup>2</sup> On Demand SoundPoint

3D Live Speech Mapping Speech Mapping Data Logging

Self Learning Swap Fit Auto Path Voice Indicators

Music Memory Memories Television Memory MT Connect™

Live Real Ear Measurement Self Check T<sup>2</sup> On Demand SoundPoint

Speech Mapping Data Logging

Self Learning Swap Fit Auto Path Voice Indicators

Music Memory Memories Television Memory MT Connect™

Live Real Ear Measurement Self Check T<sup>2</sup> On Demand SoundPoint

Speech Mapping Data Logging

# mobility™

VIBRANT 50

Designed for vibrant lifestyles, includes 16 channels and 16 bands for optimal high-resolution sound imaging, providing optimum flexibility and performance in a broad range of demanding listening environments that include high levels of background noise.

Premium ear-to-ear processing with Binaural Spatial Mapping.

Synchronized User Adjustments for volume and memory. Also includes Binaural Telephone Mode.

Compatibility with all SurfLink Accessories.

Strategic Feedback Control – leading technology that ensures wide fitting ranges and a fast response to more complex feedback.

Premium SpeechTrax optimizes 16-channel high resolution for maximum performance in noise.

Dynamic Directional Detection – designed to better perform in highly complex environments with background noise.

Live Real Ear Measurement measures the hearing aid output in your patient's ear, ensuring an accurate fit to prescriptive targets.

Live Speech Mapping verifies hearing aid's processing of speech, or any live acoustic input, in real time; 3D display engages patient and family in fitting process.

Intuitive Features include Music Processing, Voice Indicators, Automatic Telephone Solutions and Television Processing.

# mobility™

ACTIVE 40

Designed for active lifestyles, includes 12 channels and 12 bands for high-resolution sound imaging, providing excellent performance in a broad range of demanding listening environments that include moderate levels of background noise.

Advanced ear-to-ear processing with Binaural Spatial Mapping.

Synchronized User Adjustments for volume and memory. Also includes Binaural Telephone Mode.

Compatibility with all SurfLink Accessories.

Strategic Feedback Control – leading technology that ensures wide fitting ranges and a fast response to more complex feedback.

Advanced SpeechTrax utilizes 12-channel high resolution for successful performance in noise.

Dynamic Directional Detection – designed to better perform in highly complex environments with background noise.

Live Real Ear Measurement measures the hearing aid output in your patient's ear, ensuring an accurate fit to prescriptive targets.

Live Speech Mapping verifies hearing aid's processing of speech, or any live acoustic input, in real time.

Intuitive Features include Automatic Telephone Solutions and Television Processing.



# mobility™

SOCIAL 30

Designed for social lifestyles, includes 8 channels and 8 bands for precise sound imaging, providing good performance in environments with minimal to moderate levels of background noise.

Synchronized User Adjustments for volume and memory. Also includes Binaural Telephone Mode.

Compatibility with all SurfLink Accessories.

Strategic Feedback Control – leading technology that ensures wide fitting ranges and a fast response to more complex feedback.

Select SpeechTrax utilizes 8-channel resolution to improve performance in noise.

Dynamic Directional Detection – designed to better perform in highly complex environments with background noise.

Live Real Ear Measurement measures the hearing aid output in your patient's ear.

Live Speech Mapping verifies hearing aid's processing of speech, or any live acoustic input, in real time.

Intuitive Features include Automatic Telephone Solutions and Television Processing.